REFERENCE P MULTICALIBRATOR

5 x 1 ml

REF CP 0005 CH

LOT -----

componente component composant componente	valore di calibrazione / calibration value / valeur de calibration / valor de calibración					unità unit unité unidad
	LEVEL 1	LEVEL 2	LEVEL 3	LEVEL 4	LEVEL 5	
Albumina Albumin Albumine Albúmina	811 8.11	1622 16.22	3243 32.43	6487 64.87	12973 129.73	mg/dl g/l
α_1 -Antitripsina α_1 -Antitrypsin α_1 -Antitrypsine α_1 -Antitripsina	29 0.29	58 0.58	117 1.17	233 2.33	466 4.66	mg/dl g/l
α_2 -Macroglobulina α_2 -Macroglobulin α_2 -Macroglobuline α_2 -Macroglobulina	0.4 0.44	89 0.89	177 1.77	354 3.54	708 7.08	mg/dl g/l
α_1 -Glicoproteina acida α_1 -Acid Glycoprotein α_1 -Glycoprotéine acide α_1 -Glicoproteína ácida	16 0.16	31 0.31	63 0.63	125 1.25	250 2.50	mg/dl g/l
Ceruloplasmina Ceruloplasmin Céroliplasmine Ceruloplasmina	6 60	12 120	24 240	49 490	97 970	mg/dl mg/l
Aptoglobina Haptoglobin Haptoglobine Haptoglobina	24 0.24	48 0.48	95 0.95	191 1.91	381 3.81	mg/dl g/l
Transferrina Transferrin Transferrine Transferrina	55 0.55	111 1.11	222 2.22	443 4.43	886 8.86	mg/dl g/l
IgA	41 0.41	82 0.82	164 1.64	329 3.29	657 6.57	mg/dl g/l
IgG	200 2.00	400 4.00	800 8.00	1600 16.00	3200 32.00	mg/dl g/l
IgM	27 0.27	53 0.53	106 1.06	213 2.13	425 4.25	mg/dl g/l
Complemento C3 Complement C3 Complément C3 Complemento C3	21 0.21	41 0.41	82 0.82	164 1.64	328 3.28	mg/dl g/l
Complemento C4 Complement C4 Complément C4 Complemento C4	4 0.04	7 0.07	14 0.14	28 0.28	56 0.56	mg/dl g/l
Antitrombina III Antithrombin III Antithrombine III Antitrombina III	4 40	9 90	18 180	35 350	70 700	mg/dl mg/l

REFERENCE P MULTICALIBRATOR

5 x 1 ml

REF CP 0005 CH

LOT -----

componente component composant componente	valore di calibrazione / calibration value / valeur de calibration / valor de calibración					unità unit unité unidad
	LEVEL 1	LEVEL 2	LEVEL 3	LEVEL 4	LEVEL 5	
Inibitore Esterasi C1 C1 Esterase Inhibitor Inhibiteur Estérase C1 Inhibidor de Esterasa C1	4 40	8 80	16 160	32 320	63 630	mg/dl mg/l
Catena Leggera Kappa Kappa Light Chain Chaîne Légère Kappa Cadena Ligera Kappa	48 0.48	95 0.95	190 1.90	380 3.80	760 7.60	mg/dl g/l
Catena Leggera Lambda Lambda Light Chain Chaîne Légère Lambda Cadena Ligera Lambda	25 0.25	49 0.49	99 0.99	177 1.77	394 3.94	mg/dl g/l
Prealbumina Prealbumin Préalbumine Preiblúmina	5 50	10 100	20 200	40 400	80 800	mg/dl mg/l

Italiano

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Uso previsto

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Bibliografia

"Certification of proteins in the human serum - Certified Reference Material ERM-DA470k/IFCC", European Community report, 2008

English

HUMAN BASED CALIBRATOR SET

Intended use

Reference P Multicalibrator is composed of a set of sera, aimed at the preparation of calibration curves for quantitative determination of proteins in human serum, through immunochemical techniques.

Package components

Calibrators (liquid) 5 x 1 ml

Composition

Reference P Multicalibrator is a liquid, defibrinated and delipidated human plasma, containing stabilizers, sodium azide (< 0,1%) as preservative and the following human proteins: albumin, α1-antitrypsin, α2-macroglobulin, α1-acid glycoprotein, ceruloplasmin, haptoglobin, transferrin, IgA, IgG, IgM, complement C3, complement C4, antithrombin III, C1 esterase inhibitor, kappa light chain, lambda light chain, prealbumin.

Assigned values

Protein values have been assigned to ERM-DA470k/IFCC, a Reference Preparation for proteins in Human Serum (RPPHS). RPPHS is a certified reference material developed by the College of American Pathologists (CAP), the Bureau Communautaire de Reference and the International Federation of Clinical Chemists (IFCC).

Storage and stability

Store at 2-8°C.

The calibrators are ready for use; after vial opening, sera are stable for 6 weeks if stored tightly closed at 2-8°C. Do not freeze.

Test procedure

Use supplied materials as provided.

Additional materials required: automatic micropipette, general laboratory equipment.

Precautions and warnings

For in vitro diagnostic use only.

Exercise the normal precautions required for handling all laboratory reagents. This product has been prepared exclusively from the blood of donors tested individually and shown by FDA approved methods to be free from HBsAg and antibodies to HCV and HIV. However, because no test method can offer complete assurance that products derived from human sources will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens.

Disposal

This product is made to be used in professional laboratories. Please consult local regulations for a correct waste disposal.

P501: Dispose of contents according to national/international regulations.

References

"Certification of proteins in the human serum - Certified Reference Material ERM-DA470k/IFCC", European Community report, 2008

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